

Laboratory
Exclusive

Laboratory Product



Cell Culture Flask



Petri Dish



Inoculation Loop



Cell Spreader

With "caring for human health" as the initial intention.
With the mission of "improving the quality of life".

永岳医疗科技（昆山）有限公司

YONG YUE MEDICAL TECHNOLOGY (KUNSHAN)CO.,LTD

COMPANY PROFILE

Yong Yue Medical Technology (Kunshan) Co., Ltd. is a manufacturing enterprise specializing in research and development, production and sales of medical devices and laboratory consumables.

It has passed ISO13485 quality management system certification and European union certification 0123.

Products are popular in more than 80 countries and regions.

YONG YUE MEDICAL TECHNOLOGY (KUNSHAN)CO.,LTD

永岳医疗科技昆山有限公司





YONGYUE
永 | 岳 | 医 | 疗

Laboratory
Exclusive



Cell Culture Flask

Cell Culture Flask

Cell Culture Flask



non-cytotoxic



Non-pyrogenic



DNase & RNase free , human DNA free

Medical
grade
polystyrene

Gamma
radiation
sterilized

Available sizes include
25cm², 75cm², or 175cm²

Cell Culture Flask



25cm²



75cm²



175cm²

Name	Specs	Number	QTY
Cell Culture Flask	T25cm ² cell culture flask,vent cap	YCCF-25-VC	12 pcs/bag, 25 pcs/carton, 300 pcs/carton
	T75cm ² cell culture flask,vent cap	YCCF-75-VC	5pcs/bag, 18pcs/carton, 90pcs/carton
	T175cm ² cell culture flask,vent cap	YCCF-175-VC	5 pcs/bag, 10 pcs/box, 50 pcs/box



25cm²



75cm²



175cm²

Name	Specs	Number	QTY
Cell Culture Flask	T25cm ² cell culture flask, seal cap	YCCF-25-SC	12 pcs/bag, 25 pcs/carton, 300 pcs/carton
	T75cm ² cell culture flask, seal cap	YCCF-75-SC	5pcs/bag, 18pcs/carton, 90pcs/carton
	T175cm ² cell culture flask, seal cap	YCCF-175-SC	5 pcs/bag, 10 pcs/box, 50 pcs/box

Medical
grade
polystyrene(PP)

sterilized by
ethylene
oxide

Petri Dish



non-cytotoxic



Non-pyrogenic



DNase & RNase free , human DNA free

Various product models,available for
different demand of customers

Petri Dish

Round Culture Dishes



Special design makes it easier to stack and store, and saves storage space; Designed with flat and transparent surface, thickened materials;



30*15mm



60*15mm



70*15mm



90*15mm

Name	Specs	Number	QTY	DIM
Round Culture Dishes	35*15mm	YCD-01	10 pcs/bag, 2000 pcs/carton	42*42*27cm
	60*15mm	YCD-02	10 pcs/bag, 1000 pcs/carton	63*33*32cm
	70*15mm	YCD-03	10 pcs/bag, 1000 pcs/carton	74*38*35cm
	90*15mm	YCD-04	10 pcs/bag, 500 pcs/carton	48*48*32cm
	Automation 90*15mm	YCD-09	10 pcs/bag, 500 pcs/carton	48*48*32cm

Petri Dish

Round Culture Dishes



Special design makes it easier to stack and store, and saves storage space; Designed with flat and transparent surface, thickened materials;



90*20mm



100*15mm



120*15mm



150*15mm

Name	Specs	Number	QTY	DIM
Round Culture Dishes	90*20mm	YCD-05	10 pcs/bag, 500 pcs/carton	48*48*43cm
	100*15mm	YCD-10	10 pcs/bag, 500 pcs/carton	53*53*32cm
	120*15mm	YCD-11	10 pcs/pack, 320 pcs/carton	49*49*44cm
	150*15mm	YCD-12	10 pcs/bag, 200 pcs/carton	77*32*35cm


Petri Dish

Round Culture Dishes

Special design makes it easier to stack and store, and saves storage space; Designed with flat and transparent surface, thickened materials; Ventilation point design, convenient for gas exchange during cultivation;



2rooms 90*15mm



3rooms 90*15mm



4rooms 90*15mm

Name	Specs	Number	QTY	DIM
Round Culture Dishes	2rooms 90*15mm	YCD-06	10 pcs/bag, 500 pcs/carton	48*48*32cm
	3rooms 90*15mm	YCD-08	10 pcs/bag, 500 pcs/carton	48*48*32cm
	4rooms 90*15mm	YCD-07	10 pcs/bag, 500 pcs/carton	48*48*32cm

Petri Dish

Square Culture Dishes

Special design makes it easier to stack and store, and saves storage space; Designed with flat and transparent surface, thickened materials; Ventilation point design, convenient for gas exchange during cultivation;



100*100mm
with Non Scale



100*100mm
with Scale



130*130mm

Name	Specs	Number	QTY	DIM
Square Culture Dishes	100*100mm with Non Scale	YCD-13-NS	10 pcs/pack, 500 pcs/carton	56*56*37cm
	100*100mm with Scale	YCD-13-S	10 pcs/pack, 500 pcs/carton	56*56*37cm
	130*130mm	YCD-14	10 pcs/bag, 300 pcs/carton	67*41*35cm

Laboratory
Exclusive

Cell Culture Plate

Cell Culture Plate
Cell Culture Plate



non-cytotoxic



Non-pyrogenic



DNase & RNase free , human DNA free

Medical
grade
polystyrene(PP)

Various product models,available for
different demand of customers

Cell Culture Plate

Flat Bottom



6 well



12 well



24 well

Name	Specs	Number	QTY	DIM
Cell Culture Plate	6 well Culture Plate Flat Bottom	YCP-6-FB	200pcs/carton	43*29*43cm
	12 well Culture Plate Flat Bottom	YCP-12-FB	200pcs/carton	43*29*43cm
	24 well Culture Plate Flat Bottom	YCP-24-FB	200pcs/carton	43*29*43cm

Flat Bottom

U-Bottom


V-Bottom



96 well



96 well



96 well

Name	Specs	Number	QTY	DIM
Cell Culture Plate	96 well Culture Plate Flat Bottom	YCP-96-FB	Single pack 250/carton	45*40*30cm
	96 well Culture Plate U-Bottom	YCP-96-UB	Single pack 250/carton	45*40*30cm
	96 well Culture Plate V-Bottom	YCP-96-VB	Single pack 250/carton	45*40*30cm

Medical
grade
polystyrene(PP)

Gamma
radiation
sterilized



Inoculation Loop inoculation loop inoculation loop



non-cytotoxic



Non-pyrogenic



DNase & RNase free , human DNA free

Loops available
in 1 μ L and 10 μ L sizes
with strict measurement

Inoculation Loop

Square handle

Ultra-smooth loop surface design facilitates uniform and smooth streaking without damaging the gel surface.



Name:
Inoculation Loop
Specs: 1ul
Material: AS



Name:
Inoculation Loop
Specs: 10ul
Material: AS

Name	Specs	Number	QTY	DIM
Disposable Inoculating Loop	1ul Square handle	YILN-AS1-S1	Single package 500/medium box 5000/carton	84*27.5*44.5cm
		YILN-AS1-S5	5pcs/bag 200/middle bag 10000/carton	52*36*24cm
		YILN-AS1-S10	10 pcs/pack 100/middle bag 6000/carton	37*26*30cm
		YILN-AS1-S20	20 pcs/pack 500/medium box 10000/carton	77.5*27.5*34cm
	10ul Square handle	YILN-AS10-S1	Single/Pack 500/medium box 5000/carton	84*27.5*44.5cm
		YILN-AS10-S5	5pcs/pack 200/middle bag 10000/carton	52*36*24cm
		YILN-AS10-S10	10 pcs/pack 100/middle bag 6000/carton	37*26*30cm
		YILN-AS10-S20	20 pcs/pack 500/medium box 10000/carton	77.5*27.5*34cm

Inoculation Loop

Hexagon handle

Ultra-smooth loop surface design facilitates uniform and smooth streaking without damaging the gel surface.



Name:
Inoculation Loop
Specs: 10ul
Material: AS



Name:
Inoculation Loop
Specs: 1ul
Material: AS

Name	Specs	Number	QTY	DIM
Disposable Inoculating Loop	10ul Hexagonal handle	YILN-AB10-H1	Single/Pack 250/middle bag 2500/carton	52*36*24cm
		YILN-AB10-H5	5pcs/pack 100/middle bag 5000/carton	48*33*28cm
		YILN-AB10-H10	10 pcs/pack 100/medium box 10000/carton	84*44*29cm
		YILN-AB10-H20	20 pcs/pack 500/middle bag 5000/carton	52*36*22cm
	1ul Hexagonal handle	YILN-AB1-H1	Single pack 250/middle bag 2500/carton	52*36*24cm
		YILN-AB1-H5	5pcs/pack 100/middle bag 5000/carton	48*33*28cm
		YILN-AB1-H10	10 pcs/pack 100/medium box 10000/carton	84*44*29cm
		YILN-AB1-H20	20 pcs/pack 500/middle bag 5000/carton	52*36*22cm

Inoculation Loop

Round handle

Ultra-smooth loop surface design facilitates uniform and smooth streaking without damaging the gel surface.




Name:
Inoculation Loop
Specs: 10ul
Material: AS



Name:
Inoculation Loop
Specs: 1ul
Material: AS

Name	Specs	Number	QTY	DIM
Disposable Inoculating Loop	10ul Round handle	YILN-AB10-R1	Single pack 100/medium box 5000/carton	68*48.5*38.5cm
		YILN-AB10-R5	5pcs/pack 200/middle bag 10000/carton	52*35*30cm
		YILN-AB10-R10	10 pcs/pack 100/middle bag 10000/carton	55*35*21cm
		YILN-AB10-R20	20 pcs/pack 100/middle bag 10000/carton	55*35*25cm
	1ul Round handle	YILN-AB1-R1	Single pack 100/medium box 5000/carton	68*48.5*38.5cm
		YILN-AB1-R5	5pcs/pack 200/middle bag 10000/carton	52*35*30cm
		YILN-AB1-R10	10 pcs/pack 100/middle bag 10000/carton	55*35*21cm
		YILN-AB1-R20	20 pcs/pack 100/middle bag 10000/carton	55*35*25cm



Laboratory
Exclusive

Cell Spreader

Cell Spreader

Cell Spreader



non-cytotoxic



Non-pyrogenic



DNase & RNase free , human DNA free

Cell Spreader

L shaped Spreaders



The surface is smooth, and the tail-up design reduces damage to the agar surface; the design is easy to use.



Name:
Disposable L shaped Spreaders

Specs: 10ul
Material: PP


Name	Material	Number	QTY	DIM
Disposable L shaped Spreaders	PP	YDS-L-P1	Single pack 100/medium box 5000/carton	68*48.5*38.5cm
		YDS-L-P5	5pcs/pack 200/middle bag 10000/carton	52*35*30cm
		YDS-L-P10	10 pcs/pack 100/middle bag 10000/carton	55*35*21cm



Cell Spreader

T shaped Spreaders


The surface is smooth, and the tail-up design reduces damage to the agar surface; the design is easy to use.



Name:
Disposable L shaped Spreaders

Specs: 10ul
Material: PP

Name	Material	Number	QTY	DIM
Disposable T shaped Spreaders	PP	YDS-T-P1	Single pack 200/bag, 2000/carton	49*34*33cm
		YDS-T-P5	Package of 5 500/bag, 5000/carton	50*37*34cm
		YDS-T-P10	Pack of 10 500/middle bag, 5000/carton	51*33*36cm



Laboratory
Exclusive

Transfer Pipette

Transfer Pipette
Transfer Pipette



non-cytotoxic



Non-pyrogenic




DNase & RNase free , human DNA free

Medical
grade
polystyrene(PE)

Various product models,available for
different demand of customers

Transfer Pipette



Name	Specs	Number	QTY
Disposable Transfer Pipette	25ul,75mm	YBSP-025-S75	
	40ul,70mm	YBSP-040-S70	
	100UL,70mm	YBSP-0100-S70	Single pack 500/middle bag,5000pcs/CS
	200UL,70mm	YCCF-25-SC	20,000/box in bulk
	400UL,70mm	YCCF-75-SC	Single Pack 500/medium pack 5000/case
	500ul ,70mm	YCCF-175-SC	20,000/box in bulk

Transfer Pipette



Name	Specs	Number	QTY
Disposable Transfer Pipette	1ml,145mm	YBSP-1-S145	Single Pack 500/medium pack 2000/case
		YBSP-1-B145	Bulk 500/box 5000/box
	2ml,150mm	BS210110-0201	Single Pack 500/medium pack 2000/case
		BS210110-0202	Bulk 500/box 5000/box
	2ml,160mm	YBSP-3-S160	Single Pack 500/medium pack 2000/case
		YBSP-3-B160	Bulk 500/box 5000/box
	5ml,75mm	YBSP-5-S75	Single Pack, 100/bag, 3000pcs/case
		YBSP-5-B75	100 bags in bulk, 5,000 pcs/carton
	10ml,100mm	YBSP-10-S100	Single Pack 500/medium pack 5000/case



Cell Culture
Flask



TEL: 400-000-9961

FAX: 0512-57260163

Email: sales@yonyue.com

ADD: No.6 Jindie Road, Kunshan City,
Kunshan 215301, PR China

永岳医疗科技（昆山）有限公司

YONG YUE MEDICAL TECHNOLOGY (KUNSHAN) CO., LTD